

AP721914-05 Zaroa melegítő tasak - Használati utasítás (HUN)

Hogyan használjuk:

Hajlítsa meg a belső fémforgácsot, amíg a folyadék el nem kezd kristályosodni és felmelegedni. Óvatosan nyomja meg, hogy a kristályokat kiegyenlítse, és helyezze a kívánt helyre.

Az újrafelhasználás módja:

Merítse a zacskót teljesen forrásban lévő vízbe, amíg vissza nem cseppfolyósodik (Ez körülbelül 5 percet vesz igénybe)

Vigyázat:

Hűtse le teljesen, vagy tegye hideg vízbe 5 percre, mielőtt újra felhasználná.

Ne szúrja ki a zacskót

Ne használja elektromos vízfóralóban vagy mikrohullámú sütőben.

Cikkszám: AP721914-05

Forgalmazó: Cool collection, Budapest Pf.:176,
1476, Magyarország

Cégjegyzékszám: 01-09-276783

Származási ország: Kína

AP721914-05 Zaroa heat pack - User manual (ENG)

How to use:

Bend the inside metal chip until the liquid starts to crystallize and heat up.

Press gently to even at the crystals and place in the desired area.

How to reuse:

Fully immerse the bag into boiling water until it returns to liquid state

(This will take about 5 minutes)

Caution:

Cool completely or put into cold water for 5 minutes before reusing

Do not pierce the bag

Do not use in electric kettle or microwave.

Itemnumber: AP721914-05

Distributor: Cool collection, Budapest POB.176,
1476, Hungary

Registration Nr.: 01 09 276783

Made in China

ITEM 6783

Instrucciones de uso:

Doblar la pieza metálica hasta que el líquido comience a calentarse y cristalizarse.

Mezclar suavemente y colocar en la zona deseada.

Reutilización:

Sumergir la bolsa completamente en agua hirviendo, hasta que vuelva a estado líquido.
(Esto tardará unos 5 minutos aproximadamente)

Precaución:

Dejar enfriar completamente o meter en agua fría durante cinco minutos antes de reutilizar
No perforar la bolsa.
No utilizar calentadores eléctricos ni microondas para calentar la bolsa.

ITEM 6783

How to use:

Bend the inside metal chip until the liquid starts to crystallize and heat up.

Press gently to even at the crystals and place in the desired area.

How to reuse:

Fully immerse the bag into boiling water until it returns to liquid state
(This will take about 5 minutes)

Caution:

Cool completely or put into cold water for 5 minutes before reusing
Do not pierce the bag
Do not use in electric kettle or microwave.